


Product datasheet for **GA212929**

Mouse Taf3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Taf3 CRISPRa kit - CRISPR gene activation of mouse TATA-box binding protein associated factor 3
Symbol:	Taf3
Components:	GA212929G1, Taf3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212929G2, Taf3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212929G3, Taf3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5HZG4</u>
Synonyms:	140kDa; 4933439M23Rik; AW539625; mTAFII140; TAF140; TAFII-140; TAFII140
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)