


Product datasheet for **GA212939**

Mouse Tmem164 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem164 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 164
Symbol:	Tmem164
Components:	GA212939G1, Tmem164 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212939G2, Tmem164 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212939G3, Tmem164 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PHN7</u>
Synonyms:	AI316850; AI852450; AW547186; F730011B02
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)