

Product datasheet for **GA212943**

Mouse Tyw3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tyw3 CRISPRa kit - CRISPR gene activation of mouse tRNA-yW synthesizing protein 3 homolog (<i>S. cerevisiae</i>)
Symbol:	Tyw3
Components:	GA212943G1, Tyw3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212943G2, Tyw3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212943G3, Tyw3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BSA9
Synonyms:	5230400J09Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)