

Product datasheet for GA212985

Mouse Tbc1d31 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tbc1d31 CRISPRa kit - CRISPR gene activation of mouse TBC1 domain family, member 31
Symbol:	Tbc1d31
Components:	GA212985G1, Tbc1d31 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212985G2, Tbc1d31 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212985G3, Tbc1d31 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6NXYI</u>
Synonyms:	4B3; D330013L20Rik; Gm85; Wdr67
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

