

## Product datasheet for **GA212990**

### Mouse **Gab3** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	Gab3 CRISPRa kit - CRISPR gene activation of mouse growth factor receptor bound protein 2-associated protein 3
Symbol:	Gab3
Components:	GA212990G1, Gab3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA212990G2, Gab3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA212990G3, Gab3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q8BSM5</a></u>
Synonyms:	5930433H21Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)