

Product datasheet for **GA213007**

Mouse Trim41 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trim41 CRISPRa kit - CRISPR gene activation of mouse tripartite motif-containing 41
Symbol:	Trim41
Components:	GA213007G1, Trim41 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213007G2, Trim41 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213007G3, Trim41 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5NCC3</u>
Synonyms:	AW552703; BC020156; R75223; RINCK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)