


Product datasheet for **GA213035**

Mouse 4932415M13Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	4932415M13Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA 4932415M13 gene
Symbol:	4932415M13Rik
Components:	GA213035G1, 4932415M13Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213035G2, 4932415M13Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213035G3, 4932415M13Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	MGC130367
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)