

Product datasheet for **GA213036**

Mouse Tmem87a activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Tmem87a CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 87A |
| Symbol: | Tmem87a |
| UniProt ID: | <u>Q8BXN9</u> |
| Synonyms: | A930025J12Rik |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)