

Product datasheet for **GA213038**

Mouse Tifa activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tifa CRISPRa kit - CRISPR gene activation of mouse TRAF-interacting protein with forkhead-associated domain
Symbol:	Tifa
UniProt ID:	Q793I8
Synonyms:	T2bp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)