

Product datasheet for **GA213053**

Mouse **Ccdc73** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Ccdc73 CRISPRa kit - CRISPR gene activation of mouse coiled-coil domain containing 73
Symbol:	Ccdc73
Components:	GA213053G1, Ccdc73 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213053G2, Ccdc73 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213053G3, Ccdc73 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8CDM4</u>
Synonyms:	2210415I11Rik; 4931417A20; AV080665; BC095976
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)