

Product datasheet for **GA213057**

Mouse Tmem18 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem18 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 18
Symbol:	Tmem18
Components:	GA213057G1, Tmem18 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213057G2, Tmem18 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213057G3, Tmem18 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TUD9</u>
Synonyms:	BC024093
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)