

Product datasheet for GA213059

Mouse Syne3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Syne3 CRISPRa kit - CRISPR gene activation of mouse spectrin repeat containing, nuclear envelope family member 3
Symbol:	Syne3
UniProt ID:	Q4FZC9
Synonyms:	nesprin-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)