

Product datasheet for **GA213065**

Mouse Proser1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Proser1 CRISPRa kit - CRISPR gene activation of mouse proline and serine rich 1
Symbol:	Proser1
Components:	GA213065G1, Proser1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213065G2, Proser1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213065G3, Proser1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5PRE5</u>
Synonyms:	2810046L04Rik; 9330161F11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)