

Product datasheet for **GA213068**

Mouse Gsap activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Gsap CRISPRa kit - CRISPR gene activation of mouse gamma-secretase activating protein
Symbol:	Gsap
Components:	GA213068G1, Gsap gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213068G2, Gsap gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213068G3, Gsap gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TCV3</u>
Synonyms:	A530088I07Rik; Pion
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)