

Product datasheet for **GA213070**

Mouse Wdr25 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr25 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 25
Symbol:	Wdr25
Components:	GA213070G1, Wdr25 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213070G2, Wdr25 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213070G3, Wdr25 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	E9Q349
Synonyms:	B930090D16Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)