

Product datasheet for GA213104

Mouse Ddx46 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ddx46 CRISPRa kit - CRISPR gene activation of mouse DEAD (Asp-Glu-Ala-Asp) box polypeptide 46
Symbol:	Ddx46
UniProt ID:	Q569Z5
Synonyms:	2200005K02Rik; 8430438J23Rik; AI325430; AI957095; mKIAA0801
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)