

## Product datasheet for GA213109

### Mouse Tifab activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tifab CRISPRa kit - CRISPR gene activation of mouse TRAF-interacting protein with forkhead-associated domain, family member B
Symbol:	Tifab
UniProt ID:	<a href="#">Q8JZM6</a>
Synonyms:	MGC32345; MGC37193
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)