

Product datasheet for **GA213118**

Mouse Tnpo2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tnpo2 CRISPRa kit - CRISPR gene activation of mouse transportin 2 (importin 3, karyopherin beta 2b)
Symbol:	Tnpo2
Components:	GA213118G1, Tnpo2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213118G2, Tnpo2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213118G3, Tnpo2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q99LG2</u>
Synonyms:	1110034O24Rik; AA414969; AI464345; AI852433; IPO3; Knpb2b; Kpnb2b; TRN2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)