

Product datasheet for **GA213128**

Mouse **Wdr19** activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Wdr19 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 19 |
| Symbol: | Wdr19 |
| UniProt ID: | <u>Q3UGF1</u> |
| Synonyms: | C330027H04Rik; D330023L08Rik; DYF2; lft144; mKIAA1638; PWDMP |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)