

Product datasheet for **GA213147**

Mouse Lemd1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lemd1 CRISPRa kit - CRISPR gene activation of mouse LEM domain containing 1
Symbol:	Lemd1
Components:	GA213147G1, Lemd1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213147G2, Lemd1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213147G3, Lemd1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14C37</u>
Synonyms:	4930540I23Rik; 6430514M23Rik; A1428128; LEMP-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)