

Product datasheet for GA213148

Mouse Mylk3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mylk3 CRISPRa kit - CRISPR gene activation of mouse myosin light chain kinase 3
Symbol:	Mylk3
UniProt ID:	Q3UIZ8
Synonyms:	D830007F02Rik; MLCK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)