

## Product datasheet for **GA213150**

### Mouse P2ry10b activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	P2ry10b CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2Y, G-protein coupled 10B
Symbol:	P2ry10b
Components:	GA213150G1, P2ry10b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213150G2, P2ry10b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213150G3, P2ry10b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	AI662791; P2ry10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)