

Product datasheet for **GA213155**

Mouse Lgi3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lgi3 CRISPRa kit - CRISPR gene activation of mouse leucine-rich repeat LIG family, member 3
Symbol:	Lgi3
Components:	GA213155G1, Lgi3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213155G2, Lgi3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213155G3, Lgi3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8K406</u>
Synonyms:	AI841179; AW049851
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)