

Product datasheet for **GA213168**

Mouse Dync2li1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dync2li1 CRISPRa kit - CRISPR gene activation of mouse dynein cytoplasmic 2 light intermediate chain 1
Symbol:	Dync2li1
Components:	GA213168G1, Dync2li1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213168G2, Dync2li1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213168G3, Dync2li1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8K0T2
Synonyms:	4933404O11Rik; CGI-60; D2lic; LIC3; mD2LIC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)