

## Product datasheet for GA213170

### Mouse Slc44a3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Slc44a3 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 44, member 3
Symbol:	Slc44a3
UniProt ID:	<a href="#">Q921V7</a>
Synonyms:	BC010552
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)