

Product datasheet for **GA213177**

Mouse Plekhg1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plekhg1 CRISPRa kit - CRISPR gene activation of mouse pleckstrin homology domain containing, family G (with RhoGef domain) member 1
Symbol:	Plekhg1
Components:	GA213177G1, Plekhg1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213177G2, Plekhg1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213177G3, Plekhg1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	D10Ert733e; Gm521; mKIAA1209
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)