

Product datasheet for **GA213227**

Mouse Slc10a3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc10a3 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 10 (sodium/bile acid cotransporter family), member 3
Symbol:	Slc10a3
Components:	GA213227G1, Slc10a3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213227G2, Slc10a3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213227G3, Slc10a3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P21129</u>
Synonyms:	DXHXS253E; DXS253Eh; DXSmhG28; P3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)