

Product datasheet for **GA213228**

Mouse Tent4b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tent4b CRISPRa kit - CRISPR gene activation of mouse terminal nucleotidyltransferase 4B
Symbol:	Tent4b
Components:	GA213228G1, Tent4b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213228G2, Tent4b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213228G3, Tent4b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q68ED3</u>
Synonyms:	5730445M16Rik; 5830428A09
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)