


Product datasheet for **GA213229**

Mouse L3mbtl2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	L3mbtl2 CRISPRa kit - CRISPR gene activation of mouse L3MBTL2 polycomb repressive complex 1 subunit
Symbol:	L3mbtl2
Components:	GA213229G1, L3mbtl2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213229G2, L3mbtl2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213229G3, L3mbtl2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P59178</u>
Synonyms:	4732493N06Rik; M4mbt
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)