

Product datasheet for **GA213235**

Mouse Zfp879 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp879 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 879
Symbol:	Zfp879
Components:	GA213235G1, Zfp879 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213235G2, Zfp879 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213235G3, Zfp879 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BI99</u>
Synonyms:	Znf879
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)