

Product datasheet for GA213236

Mouse Sertad4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sertad4 CRISPRa kit - CRISPR gene activation of mouse SERTA domain containing 4
Symbol:	Sertad4
Components:	GA213236G1, Sertad4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213236G2, Sertad4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213236G3, Sertad4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A7DTG3</u>
Synonyms:	4833404A09; C130018M11Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)