


## Product datasheet for **GA213237**

### Mouse Syde2 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Syde2 CRISPRa kit - CRISPR gene activation of mouse synapse defective 1, Rho GTPase, homolog 2 (C. elegans)                                                                                                                                                                                                |
| Symbol:       | Syde2                                                                                                                                                                                                                                                                                                      |
| Components:   | GA213237G1, Syde2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA213237G2, Syde2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA213237G3, Syde2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | C430017H16Rik                                                                                                                                                                                                                                                                                              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)