

Product datasheet for **GA213238**

Mouse Zfp609 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp609 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 609
Symbol:	Zfp609
Components:	GA213238G1, Zfp609 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213238G2, Zfp609 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213238G3, Zfp609 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BZ47
Synonyms:	9830165N24; BC016271; mKIAA0295; Znf609
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)