

Product datasheet for **GA213245**

Mouse Slc39a2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc39a2 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 39 (zinc transporter), member 2
Symbol:	Slc39a2
Components:	GA213245G1, Slc39a2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213245G2, Slc39a2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213245G3, Slc39a2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	F730005G13Rik; Gm289; Gm1789; Zip2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)