


Product datasheet for **GA213255**

Mouse Prl3d3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Prl3d3 CRISPRa kit - CRISPR gene activation of mouse prolactin family 3, subfamily d, member 3
Symbol:	Prl3d3
Components:	GA213255G1, Prl3d3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213255G2, Prl3d3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213255G3, Prl3d3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Ghd6; PL-Ig; Plig
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)