

Product datasheet for **GA213271**

Mouse Slc36a3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc36a3 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 36 (proton/amino acid symporter), member 3
Symbol:	Slc36a3
Components:	GA213271G1, Slc36a3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213271G2, Slc36a3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213271G3, Slc36a3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q811P0</u>
Synonyms:	PAT3; TRAMD2; tramdorin2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)