

Product datasheet for **GA213315**

Mouse Zfp365 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp365 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 365
Symbol:	Zfp365
Components:	GA213315G1, Zfp365 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213315G2, Zfp365 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213315G3, Zfp365 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BG89</u>
Synonyms:	AI839779; AV340892; DBZ; Su48; talanin; Znf365
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)