

Product datasheet for **GA213317**

Mouse Ybey activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ybey CRISPRa kit - CRISPR gene activation of mouse ybeY metalloproteinase
Symbol:	Ybey
Components:	GA213317G1, Ybey gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213317G2, Ybey gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213317G3, Ybey gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8CAV0</u>
Synonyms:	1110068N23Rik; A130042E20Rik; C21orf57; ORF57; ORF66
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)