

Product datasheet for **GA213325**

Mouse Plppr3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plppr3 CRISPRa kit - CRISPR gene activation of mouse phospholipid phosphatase related 3
Symbol:	Plppr3
UniProt ID:	Q7TPB0
Synonyms:	Lppr3; Prg2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)