

Product datasheet for **GA213344**

Mouse Tspan8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tspan8 CRISPRa kit - CRISPR gene activation of mouse tetraspanin 8
Symbol:	Tspan8
UniProt ID:	<u>Q8R3G9</u>
Synonyms:	C76990; E330007O21Rik; Tm4sf3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)