

Product datasheet for GA213357

Mouse Pik3ip1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pik3ip1 CRISPRa kit - CRISPR gene activation of mouse phosphoinositide-3-kinase interacting protein 1
Symbol:	Pik3ip1
UniProt ID:	Q7TMJ8
Synonyms:	1500004A08Rik; 5830455E04Rik; Crkd; Hgfl
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)