


Product datasheet for **GA213358**

Mouse Ccdc157 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ccdc157 CRISPRa kit - CRISPR gene activation of mouse coiled-coil domain containing 157
Symbol:	Ccdc157
Components:	GA213358G1, Ccdc157 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213358G2, Ccdc157 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213358G3, Ccdc157 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5SPX1</u>
Synonyms:	4930562D19Rik; mKIAA1656
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)