


Product datasheet for **GA213378**

Mouse Iba57 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Iba57 CRISPRa kit - CRISPR gene activation of mouse IBA57 homolog, iron-sulfur cluster assembly
Symbol:	Iba57
Components:	GA213378G1, Iba57 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213378G2, Iba57 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213378G3, Iba57 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8CAK1</u>
Synonyms:	4930543L23Rik; A230051G13Rik; C1orf69
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)