

Product datasheet for **GA213380**

Mouse Prss38 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Prss38 CRISPRa kit - CRISPR gene activation of mouse protease, serine 38
Symbol:	Prss38
Components:	GA213380G1, Prss38 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213380G2, Prss38 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213380G3, Prss38 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3UKY7</u>
Synonyms:	Gm249; Mpn2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)