

## Product datasheet for **GA213388**

### Mouse **Arhgap44** activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Arhgap44 CRISPRa kit - CRISPR gene activation of mouse Rho GTPase activating protein 44
Symbol:	Arhgap44
UniProt ID:	<u><a href="#">Q5SSM3</a></u>
Synonyms:	6330543G20; AI840762; AU040829; AW493732; Rich2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)