

Product datasheet for **GA213389**

Mouse Usp43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Usp43 CRISPRa kit - CRISPR gene activation of mouse ubiquitin specific peptidase 43
Symbol:	Usp43
Components:	GA213389G1, Usp43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213389G2, Usp43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213389G3, Usp43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BUM9</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)