


Product datasheet for **GA213431**

Mouse Cisd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cisd3 CRISPRa kit - CRISPR gene activation of mouse CDGSH iron sulfur domain 3
Symbol:	Cisd3
Components:	GA213431G1, Cisd3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213431G2, Cisd3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213431G3, Cisd3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	B1AR13
Synonyms:	Mel-13; Mel13
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)