


Product datasheet for **GA213436**

Mouse Plekhh3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plekhh3 CRISPRa kit - CRISPR gene activation of mouse pleckstrin homology domain containing, family H (with MyTH4 domain) member 3
Symbol:	Plekhh3
Components:	GA213436G1, Plekhh3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213436G2, Plekhh3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213436G3, Plekhh3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8VCE9</u>
Synonyms:	BC020025
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)