

Product datasheet for GA213438

Mouse Tmem106a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem106a CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 106A
Symbol:	Tmem106a
UniProt ID:	<u>Q8VC04</u>
Synonyms:	0610008L10Rik; AI043106; BC022145
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)