

Product datasheet for **GA213447**

Mouse Ace3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ace3 CRISPRa kit - CRISPR gene activation of mouse angiotensin I converting enzyme (peptidyl-dipeptidase A) 3
Symbol:	Ace3
Components:	GA213447G1, Ace3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA213447G2, Ace3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA213447G3, Ace3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>D0G895</u>
Synonyms:	EG217246
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)